BRCA+16 GENES

**Hereditary predisposition to breast,
ovarian and endometrial cancer**

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Index

- **Breast and ovarian cancer**
- **BRCA⁺¹⁶ GENES**
 - Genes panel and associated risk
 - Analysis technology
 - Recommendations from medical associations
 - Indications
 - Important considerations
 - Advantages
 - Sample collection kit
 - Test Requisition Form and Informed Consent
 - Results
 - Key points
- **Contact**

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- Genes panel and associated risk

- Analysis technology

- Recommendations from medical associations

- Indications

- Important considerations

- Advantages

- Sample collection kit

- Test Requisition Form and Informed Consent

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
Breast and ovarian cancer

According to the International Agency for Research on Cancer:¹

- **Breast cancer** is the most common type of cancer in women worldwide.
- It mainly affects women, although it can also affect men.

A quarter of the tumours detected in women are breast tumours.²


*Estimated number of incident cases, females, worldwide
(top 10 cancer sites) in 2012*

1. GLOBOCAN 2008. Cancer Incidence and Mortality Worldwide: IARC Cancer Base No. 10 [Internet]. Lyon, France: International Agency for Research on Cancer; 2010.
2. Spanish Association Against Cancer (AECC).

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1. GLOBOCAN 2008. Cancer Incidence and Mortality Worldwide: IARC Cancer Base No. 10 [Internet]. Lyon, France: International Agency for Research on Cancer; 2010.
2. Spanish Association Against Cancer (AECC).

Breast and ovarian cancer

BREAST CANCER

More than **1,600,000 new cases** are diagnosed per year worldwide.¹

Approximately **1 in 8 women** will be diagnosed with breast cancer during their **lifetime**.²

5-10% of cases are **hereditary**.³

OVARIAN CANCER

Over **230,000 new cases** diagnosed per year worldwide.¹

Approximately **20%** are **hereditary**.³

1. Ferlay J, Soerjomataram I, Dikshit R, et al. Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012. *Int J Cancer* 2015;136:E359-86.
2. European commission initiative on breast cancer website: <http://ecibc.jrc.ec.europa.eu/recommendations/>
3. Nielsen FC et al. Hereditary breast and ovarian cancer: new genes in confined pathways. *Nature Reviews*. 2016;16:599-612.

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1. Ferlay J, Soerjomataram I, Dikshit R, et al. Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012. *Int J Cancer* 2015;136:E359-86.
2. European commission initiative on breast cancer website: <http://ecibc.jrc.ec.europa.eu/recommendations/>
3. Nielsen FC et al. Hereditary breast and ovarian cancer: new genes in confined pathways. *Nature Reviews*. 2016;16:599-612.

Breast and ovarian cancer: Risk factors


Lifestyle-related Risk Factors

-  Drinking alcohol (breast)
-  Being overweight or obese (both)
-  Not being physically active (breast)
-  Not having children (both)
-  Not breastfeeding (breast)
-  Birth control (breast)
-  Fertility treatment (ovarian)
-  Hormone therapy after menopause (both)

Risk Factors You Cannot Change

-  Being a woman (both)
-  Getting older (both)
-  Certain inherited genes (both)
-  Family cancer syndromes (both)
-  Family history of breast cancer (breast)
-  Personal history of breast cancer (both)
-  Early menarche before age 12 (breast)
-  Late menopause after age 55 (breast)

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
<https://www.cancer.org/cancer/>

Breast and ovarian cancer: Hereditary cancer

- **Between 5 to 10%** of cases of **breast cancer** and **20%** of **ovarian cancer** are associated with a **hereditary** nature.
- Characterised by its incidence at **early ages**, even before 40 years of age.
- Frequently associated to mutations in **BRCA1** and **BRCA2** **genes, tumour suppressor genes** involved in maintaining DNA integrity.
- There are other genes related with these cancer types that must be studied since it is estimated that only **around 25% of hereditary breast and ovarian cancer cases are due to mutations in BRCA1 and BRCA2 genes.**¹

Approximately **50% of women with mutations in the BRCA1 or BRCA2 genes do not have a history of breast or ovarian cancer.**²

1. Nielsen FC et al. Hereditary breast and ovarian cancer: new genes in confined pathways. *Nature Reviews*. 2016;16:599-612.

2. King MC, Levy-Lahad E, Lahad A. Population-Based Screening for BRCA1 and BRCA2: 2014 Lasker Award. *JAMA*. 2014;312(11):1091-2.

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
1. Nielsen FC et al. Hereditary breast and ovarian cancer: new genes in confined pathways. *Nature Reviews*. 2016;16:599-612.

2. King MC, Levy-Lahad E, Lahad A. Population-Based Screening for BRCA1 and BRCA2: 2014 Lasker Award. *JAMA*. 2014;312(11):1091-2.

Breast and ovarian cancer: Hereditary cancer


- Cumulative risk of developing breast and ovarian cancer throughout life in people with and without *BRCA1* and *BRCA2* mutations.

Modified from: Petrucelli N. et al. *BRCA1- and BRCA2-Associated Hereditary Breast and Ovarian Cancer*. GeneReviews. Last Update: December 15, 2016. Available at: www.ncbi.nlm.nih.gov/books/NBK1247/

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Modified from: Petrucelli N. et al. *BRCA1- and BRCA2-Associated Hereditary Breast and Ovarian Cancer*. GeneReviews. Last Update: December 15, 2016. Available at: www.ncbi.nlm.nih.gov/books/NBK1247/

Breast and ovarian cancer: Hereditary cancer

International experts recommend the screening of *BRCA1* and *BRCA2* genes in women \geq 30 years in the course of routine medical care.¹

VIEWPOINT Population-Based Screening for *BRCA1* and *BRCA2* 2014 Lasker Award

Mary-Claire King, PhD
Departments of
Medicine and Genome
Sciences, University of
Washington, Seattle

Ephrat Levy-Lahad,
MD
Medical Genetics
Institute, Shaare Zedek
Medical Center

The 2014 Lasker-Koshland Special Achievement Award in Medical Science has been presented to Dr Mary-Claire King to recognize and honor her "far bold and imaginative contributions to medical science and society—exemplified by her discovery of a single gene *BRCA1* that causes a... form of hereditary breast cancer..." This Viewpoint describes the application of that discovery, and suggests that population-based screening of women for *BRCA1* and *BRCA2* should become a routine part of clinical practice.

on the data then available. However, a just-completed study now provides evidence that supports offering *BRCA1* and *BRCA2* sequencing to all women.
To determine cancer risks to *BRCA1* and *BRCA2* mutation carriers identified from the general population, we conducted a study of population-based screening in the Ashkenazi Jewish population of Israel.⁶ This population was selected because its limited number of cancer-predisposing *BRCA1* and *BRCA2* mutations made the study feasible to achieve meaningful representation of all

Based on our 20 years' experience working with families with cancer-predisposing mutations in *BRCA1* and *BRCA2*, it is time to offer genetic screening of these genes to every woman

Mary-Claire King PhD; Ephrat Levy-Lahad MD; Amnon Lahad MPH

screening of these genes to every woman

to every woman, at about age 30, in the course of routine medical care. Women with cancer-predisposing mutations in *BRCA1* and *BRCA2* are a high-risk group in whom special screening and counseling can be focused.

World Health Organization criteria for population screening for genetic predisposition to disease are that the disease is an important public health burden in the target population; that the risk of disease due to mutations in the screened genes is known; and that effective interventions exist to reduce morbidity and mortality among genetically susceptible individuals.⁶ At present, the US Preventive Services Task Force (USPSTF) supports *BRCA1* and *BRCA2* testing based on family history and ancestry, but not for the entire female population, given the lack of data on risks for mutation carriers ascertained from the general population, rather than through a personal or family history of cancer.⁶ This position was correct based

on the data then available. However, a just-completed study now provides evidence that supports offering *BRCA1* and *BRCA2* sequencing to all women.
To determine cancer risks to *BRCA1* and *BRCA2* mutation carriers identified from the general population, we conducted a study of population-based screening in the Ashkenazi Jewish population of Israel.⁶ This population was selected because its limited number of cancer-predisposing *BRCA1* and *BRCA2* mutations made the study feasible to achieve meaningful representation of all

jama.com

JAMA September 17, 2014 Volume 312, Number 11 1091

Copyright 2014 American Medical Association. All rights reserved.

Downloaded From: <http://jama.jamanetwork.com/> by a Penn State Milton S Eisenhower Med Ctr User on 09/26/2014

World Health Organization criteria for population screening for genetic predisposition to disease are:

- the disease is an **important public health burden in the target population**
- that the **risk of disease** due to mutations in the screened genes **is known**
- and that **effective interventions exist to reduce morbidity and mortality** among genetically susceptible individuals

1. King MC, Levy-Lahad E, Lahad A. Population-Based Screening for *BRCA1* and *BRCA2*: 2014 Lasker Award. *JAMA*. 2014;312(11):1091-2.

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VIEWPOINT Population-Based Screening for *BRCA1* and *BRCA2*
2014 Lasker Award

Mary-Claire King, PhD
Departments of
Medicine and Genome
Sciences, University of
Washington, Seattle

**Ephrat Levy-Lahad,
MD**
Medical Genetics
Institute, Sheba Daskin
Medical Center

The 2014 Lasker-Koshland Special Achievement Award in Medical Science has been presented to Dr Mary-Claire King to recognize and honor her "far bold and imaginative contributions to medical science and society—exemplified by her discovery of a single gene *BRCA1* that causes a... form of hereditary breast cancer..." This Viewpoint describes the application of that discovery, and suggests that population-based screening of women for *BRCA1* and *BRCA2* should become a routine part of clinical practice.

on the data then available. However, a just-completed study now provides evidence that supports offering *BRCA1* and *BRCA2* sequencing to all women.⁶

To determine cancer risks to *BRCA1* and *BRCA2* mutation carriers identified from the general population, we conducted a study of population-based screening in the Ashkenazi Jewish population of Israel.⁶ This population was selected because its limited number of cancer-predisposing *BRCA1* and *BRCA2* mutations made the study feasible to achieve meaningful representation of all

Based on our 20 years' experience working with families with cancer-predisposing mutations in *BRCA1* and *BRCA2*, it is time to offer genetic screening of these genes to every woman

Mary-Claire King PhD; Ephrat Levy-Lahad MD; Amnon Lahad MPH

screening of these genes to every woman

to every woman, at about age 30, in the course of routine medical care. Women with cancer-predisposing mutations in *BRCA1* and *BRCA2* are a high-risk group in whom special screening and counseling can be focused.

World Health Organization criteria for population screening for genetic predisposition to disease are that the disease is an important public health burden in the target population; that the risk of disease due to mutations in the screened genes is known; and that effective interventions exist to reduce morbidity and mortality among genetically susceptible individuals.⁶ At present, the US Preventive Services Task Force (USPSTF) supports *BRCA1* and *BRCA2* testing based on family history and ancestry, but not for the entire female population, given the lack of data on risks for mutation carriers ascertained from the general population, rather than through a personal or family history of cancer.⁶ This position was correct based

ovarian cancer was 60% (47%) by age 60 and 83% (47%) by age 80. For *BRCA2* mutation carriers, risk was 33% (19%) by age 60 and 76% (43%) by age 80. Furthermore, these risks were significantly higher, at every age, among women born more recently than among women born earlier; a birth cohort effect also seen in prior studies. This trend likely reflects increasing prevalence of nongenetic risk factors for breast cancer, including earlier age of menarche and later ages of childbearing, factors related to improved nutrition and education for women in modern society. Notably, 50% of families found to harbor *BRCA1* or *BRCA2* mutations had no history of breast or ovarian cancer that would have triggered clinical attention. Female mutation carriers from these low-cancer-incidence families had similar cancer risks to female carriers from families with high cancer incidence. Low-cancer-incidence families were simply smaller, with fewer females who inherited *BRCA1* or *BRCA2* mutations, and hence fewer females

Corresponding
Author: Mary-Claire
King, PhD,
Departments of
Medicine and Genome
Sciences, University of
Washington, Health
Sciences, Room K-360,
Seattle, WA 98195;
7702 mcking@uw
.edu

jama.com JAMA September 17, 2014 Volume 312, Number 11 1091

Copyright 2014 American Medical Association. All rights reserved.

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Breast and ovarian cancer: Hereditary cancer

International experts recommend the screening of *BRCA1* and *BRCA2* genes in women \geq 30 years in the course of routine medical care.¹

VIEWPOINT Population-Based Screening for *BRCA1* and *BRCA2*
2014 Lasker Award

Mary-Claire King, PhD
Departments of
Medicine and Genome
Sciences, University of
Washington, Seattle

**Ephrat Levy-Lahad,
MD**
Medical Genetics
Institute, Shaare Zedek
Medical Center

The 2014 Lasker-Koshland Special Achievement Award in Medical Science has been presented to Dr Mary-Claire King to recognize and honor her "far bold and imaginative contributions to medical science and society — exemplified by her discovery of a single gene *BRCA1* that causes a ... form of hereditary breast cancer ...". This Viewpoint describes the application of that discovery, and suggests that population-based screening of women for *BRCA1* and *BRCA2* should become a routine part of clinical practice.

on the data then available. However, a just-completed study now provides evidence that supports offering *BRCA1* and *BRCA2* sequencing to all women.
To determine cancer risks to *BRCA1* and *BRCA2* mutation carriers identified from the general population, we conducted a study of population-based screening in the Ashkenazi Jewish population of Israel.⁶ This population was selected because its limited number of cancer-predisposing *BRCA1* and *BRCA2* mutations made the study feasible to achieve meaningful representation of all

Based on our 20 years' experience working with families with cancer-predisposing mutations in *BRCA1* and *BRCA2*, it is time to offer genetic screening of these genes to every woman

Mary-Claire King PhD; Ephrat Levy-Lahad MD; Amnon Lahad MPH

screening of these genes to every woman

to every woman, at about age 30, in the course of routine medical care. Women with cancer-predisposing mutations in *BRCA1* and *BRCA2* are a high-risk group in whom special screening and counseling can be focused.
World Health Organization criteria for population screening for genetic predisposition to disease are that the disease is an important public health burden in the target population; that the risk of disease due to mutations in the screened genes is known; and that effective interventions exist to reduce morbidity and mortality among genetically susceptible individuals.⁶ At present, the US Preventive Services Task Force (USPSTF) supports *BRCA1* and *BRCA2* testing based on family history and ancestry, but not for the entire female population, given the lack of data on risks for mutation carriers ascertained from the general population, rather than through a personal or family history of cancer.⁶ This position was correct based

ovarian cancer was 60% (a 7% by age 60 and 83% (a 7% by age 80. For *BRCA2* mutation carriers, risk was 33% (a 9% by age 60 and 76% (a 13% by age 80. Furthermore, these risks were significantly higher, at every age, among women born more recently than among women born earlier; a birth cohort effect also seen in prior studies. This trend likely reflects increasing prevalence of nongenetic risk factors for breast cancer, including earlier age of menarche and later ages of childbearing, factors related to improved nutrition and education for women in modern society. Notably, 50% of families found to harbor *BRCA1* or *BRCA2* mutations had no history of breast or ovarian cancer that would have triggered clinical attention. Female mutation carriers from these low-cancer-incidence families had similar cancer risks to female carriers from families with high cancer incidence. Low-cancer-incidence families were simply smaller, with fewer females who inherited *BRCA1* or *BRCA2* mutations, and hence fewer females

Corresponding
Author: Mary-Claire
King, PhD,
Departments of
Medicine and Genome
Sciences, University of
Washington, Health
Sciences, Room K-360,
Seattle, WA 98195;
7702 mckingaw
@u.wa.edu

JAMA September 17, 2014 Volume 312, Number 11 1091

Copyright 2014 American Medical Association. All rights reserved.

- **Low-cancer-incidence families were simply smaller, with fewer females who inherited *BRCA1* or *BRCA2* mutations, and hence fewer females who developed breast or ovarian cancer. Absent population-wide screening, women with *BRCA1* or *BRCA2* mutations from such families would not have been identified until they developed cancer.**
- As population-based screening for *BRCA1* and *BRCA2* among adult women becomes a routine part of clinical practice, **other genes are expected to be phased into the process.**

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VIEWPOINT Population-Based Screening for *BRCA1* and *BRCA2*
2014 Lasker Award

Mary-Claire King, PhD
Departments of
Medicine and Genome
Sciences, University of
Washington, Seattle

**Ephrat Levy-Lahad,
MD**
Medical Genetics
Institute, Shaare Zedek
Medical Center

The 2014 Lasker-Koshland Special Achievement Award in Medical Science has been presented to Dr Mary-Claire King to recognize and honor her "far bold and imaginative contributions to medical science and society... exemplified by her discovery of a single gene *BRCA1* that causes a... form of hereditary breast cancer..." This Viewpoint describes the application of that discovery, and suggests that population-based screening of women for *BRCA1* and *BRCA2* should become a routine part of clinical practice.

on the data then available. However, a just-completed study now provides evidence that supports offering *BRCA1* and *BRCA2* sequencing to all women.
To determine cancer risks to *BRCA1* and *BRCA2* mutation carriers identified from the general population, we conducted a study of population-based screening in the Ashkenazi Jewish population of Israel.⁶ This population was selected because its limited number of cancer-predisposing *BRCA1* and *BRCA2* mutations made the study feasible to achieve meaningful representation of all

Based on our 20 years' experience working with families with cancer-predisposing mutations in *BRCA1* and *BRCA2*, it is time to offer genetic screening of these genes to every woman

Mary-Claire King PhD; Ephrat Levy-Lahad MD; Amnon Lahad MPH

screening of these genes to every woman

to every woman, at about age 30, in the course of routine medical care. Women with cancer-predisposing mutations in *BRCA1* and *BRCA2* are a high-risk group in whom special screening and counseling can be focused.
World Health Organization criteria for population screening for genetic predisposition to disease are that the disease is an important public health burden in the target population; that the risk of disease due to mutations in the screened genes is known; and that effective interventions exist to reduce morbidity and mortality among genetically susceptible individuals.⁶ At present, the US Preventive Services Task Force (USPSTF) supports *BRCA1* and *BRCA2* testing based on family history and ancestry, but not for the entire female population, given the lack of data on risks for mutation carriers ascertained from the general population, rather than through a personal or family history of cancer.⁶ This position was correct based

ovarian cancer was 60% (47%) by age 60 and 83% (47%) by age 80. For *BRCA2* mutation carriers, risk was 33% (19%) by age 60 and 76% (43%) by age 80. Furthermore, these risks were significantly higher, at every age, among women born more recently than among women born earlier; a birth cohort effect also seen in prior studies. This trend likely reflects increasing prevalence of nongenetic risk factors for breast cancer, including earlier age of menarche and later ages of childbearing, factors related to improved nutrition and education for women in modern society. Notably, 50% of families found to harbor *BRCA1* or *BRCA2* mutations had no history of breast or ovarian cancer that would have triggered clinical attention. Female mutation carriers from these low-cancer-incidence families had similar cancer risks to female carriers from families with high cancer incidence. Low-cancer-incidence families were simply smaller, with fewer females who inherited *BRCA1* or *BRCA2* mutations, and hence fewer females



Corresponding
Author: Mary-Claire
King, PhD,
Departments of
Medicine and Genome
Sciences, University of
Washington, Health
Sciences, Room 4-360,
Seattle, WA 98195;
7702 mcking@uw
.edu

jama.com JAMA September 17, 2014 Volume 312, Number 11 1091


Copyright 2014 American Medical Association. All rights reserved.

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Breast and ovarian cancer: Hereditary cancer

External factors are especially important in people carrying mutations since they may alter the unaffected gene

- Every individual inherits 2 copies of each gene, one from the mother and the other from the father.
- In BRCA, if we inherit altered copy, only the not affected gene will protect the cell from tumour development (heterocygous).
- People with only one altered BRCA copy have a greater predisposition to tumour development.



Facing the exposure to different risk factors, the no mutated gene is altered and lose its functionality

In these situations, the probability of developing a tumour is very high.

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Genes panel and associated risk

Analysis technology

Recommendations from medical associations

Indications

Important considerations

Advantages

Sample collection kit

Test Requisition Form and Informed Consent

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- Genes **BRCA1, BRCA2 + 16 genes** associated to breast, ovarian and endometrial cancer.
- **Designed and developed by genetic SYNLAB group experts, in line with the NCCN guidelines**, including the most relevant genes for which the guidelines recommend a specific patient management.
- **NGS sequencing** (*Next Generation Sequencing*) **with paired-end reads** of the genes of the panel, which allows to detect any pathogenic mutation or variant of uncertain significance (VUS).
- **Large deletions and duplications analysis** in *BRCA1, BRCA2* and *EPCAM* genes through MLPA (*Multiplex Ligation Probe Amplification*).
- **Pathogenic and probably pathogenic mutations** are confirmed by **Sanger sequencing**.

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BRCA+16GENES

- The genes included in **BRCA+16 GENES** are involved in **cell cycle control and DNA repair** during cell division.
- Mutations in these genes lead to a **loss of cell control and capacity for DNA repair**, which may implies a greater risk of developing cancer than the general population.
- Abnormalities in the genes included in the **BRCA+16 GENES** panel represent an **increased risk of suffering from hereditary breast, ovarian and endometrial cancer**.

ANALYTIC PERFORMANCE

SENSITIVITY	SPECIFICITY	MINIMUM COVERAGE
99%	>96%	20x

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BRCA+16GENES: Gene panel and associated risk

HIGH RISK GENES	LIFETIME RISK	OTHER CANCERS
<i>BRCA1</i>	38-84% breast cancer	Ovarian, prostate, pancreas
<i>BRCA2</i>	38-84% breast cancer	Ovarian, prostate, pancreas
<i>TP53</i> (Li Fraumeni syndrome)	≤79% breast cancer	Gastric, sarcoma, brain tumor
<i>PTEN</i> (Cowden syndrome)	25-50% breast cancer	Melanoma, prostate, endometrium
<i>CDH1</i> (Hereditary diffuse gastric cancer)	39-52% breast cancer	Gastric
<i>STK11</i> (Peutz-Jeghers syndrome)	32-54% breast cancer	Pancreas, gastrointestinal, sex cord–gonadal stromal tumour

MODERATE RISK GENES	LIFETIME RISK	OTHER CANCERS
<i>PALB2</i>	44% breast cancer	Pancreas
<i>CHEK2</i>	32% breast cancer	Pancreas, lung
<i>ATM</i>	30% breast cancer	Pancreas
<i>NBN</i>	30% breast cancer	Ovarian

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BRCA+16GENES: Gene panel and associated risk

REPAIR GENES	LIFETIME RISK	OTHER CANCERS
<i>RAD51C</i>	6,1% ovarian cancer	Breast
<i>RAD51D</i>	13,5% ovarian cancer	Breast
<i>BRIP1</i>	4-12% ovarian cancer	Breast

LYNCH SYNDROME	LIFETIME RISK	OTHER CANCERS
<i>MLH1, MSH2, MSH6, PMS2, EPCAM</i>	4-12% ovarian cancer	Colorectal
	16-60% endometrial cancer	

REFERENCES

1. Easton DF, Pharoah PD, Antoniou AC, et al. Gene-panel sequencing and the prediction of breast-cancer risk. *N Engl J Med.* 2015;372(23):2243-57
2. Kurian AW, Hare EE, Mills MA, et al. Clinical evaluation of a multiple-gene sequencing panel for hereditary cancer risk assessment. *J Clin Oncol.* 2014;32(19):2001-9.
3. Tung N, Domchek SM, Stadler Z. Counselling framework for moderate-penetrance cancer-susceptibility mutations. *Nat Rev Clin Oncol.* 2016;13(9):581-8.
4. Giardiello FM, Brensinger JD, Tersmette AC, et al. Very High Risk of Cancer in Familial Peutz-Jeghers Syndrome. *GASTROENTEROLOGY* 2000;119:1447–1453.
5. Møller P, Seppälä TT, Bernstein I, et al. Cancer risk and survival in path_MMR carriers by gene and gender up to 75 years of age: a report from the Prospective Lynch Syndrome Database. *Gut.* 2017. doi: 10.1136/gutjnl-2017-314057.
6. Petrucelli N. et al. BRCA1- and BRCA2-Associated Hereditary Breast and Ovarian Cancer. *GeneReviews.* Last Update: December 15, 2016. Available at: www.ncbi.nlm.nih.gov/books/NBK1247
7. Provenzale D, Gupta S, Ahnen DJ, et al. Genetic/Familial High-Risk Assessment: Colorectal Version 1.2016, NCCN Clinical Practice Guidelines in Oncology. *J Natl Compr Canc Netw.* 2016;14(8):1010-30.
8. OMIM databases: Online Mendelian Inheritance in Man@ An Online Catalog of Human Genes and Genetic Disorders - <https://www.omim.org/>

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
REFERENCES

1. Easton DF, Pharoah PD, Antoniou AC, et al. Gene-panel sequencing and the prediction of breast-cancer risk. *N Engl J Med.* 2015;372(23):2243-57
2. Kurian AW, Hare EE, Mills MA, et al. Clinical evaluation of a multiple-gene sequencing panel for hereditary cancer risk assessment. *J Clin Oncol.* 2014;32(19):2001-9.
3. Tung N, Domchek SM, Stadler Z. Counselling framework for moderate-penetrance cancer-susceptibility mutations. *Nat Rev Clin Oncol.* 2016;13(9):581-8.
4. Giardiello FM, Brensinger JD, Tersmette AC, et al. Very High Risk of Cancer in Familial Peutz-Jeghers Syndrome. *GASTROENTEROLOGY* 2000;119:1447–1453.
5. Møller P, Seppälä TT, Bernstein I, et al. Cancer risk and survival in path_MMR carriers by gene and gender up to 75 years of age: a report from the Prospective Lynch Syndrome Database. *Gut.* 2017. doi: 10.1136/gutjnl-2017-314057.
6. Petrucelli N. et al. BRCA1- and BRCA2-Associated Hereditary Breast and Ovarian Cancer. *GeneReviews.* Last Update: December 15, 2016. Available at: www.ncbi.nlm.nih.gov/books/NBK1247
7. Provenzale D, Gupta S, Ahnen DJ, et al. Genetic/Familial High-Risk Assessment: Colorectal Version 1.2016, *NCCN Clinical Practice Guidelines in Oncology. J Natl Compr Canc Netw.* 2016;14(8):1010-30.
8. OMIM databases: Online Mendelian Inheritance in Man@ An Online Catalog of Human Genes and Genetic Disorders - <https://www.omim.org/>

BRCA+16GENES: International Guidelines

9

NCCN Guidelines® Insights Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

Featured Updates to the NCCN Guidelines

Mary B. Daly, MD, PhD^{1,2}; Robert Pilarski, MS, CGC³; Michael Berry, MD⁴; Sandra S. Buys, MD⁵; Meagan Farmer, MS, CGC⁶; Susan Friedman, DVM⁷; Judy E. Garber, MD, MPH⁸; Noah D. Kauff, MD⁹; Seema Khan, MD¹⁰; Catherine Klein, MD¹¹; Wendy Kohlmann, MS, CGC¹²; Allison Kurian, MD, MS¹³; Jennifer K. Litton, MD¹⁴; Lisa Madlensky, PhD, CGC¹⁵; Sofia D. Merajver, MD, PhD¹⁶; Kenneth Offit, MD¹⁷; Toyo Pal, MD¹⁸; Gwen Reiser, MS, CGC¹⁹; Kristen Mahoney Shannon, MS, CGC²⁰; Elizabeth Swisher, MD²¹; Shweta Vinayak, MD²²; Nicoleta C. Voian, MD, MPH²³; Jeffrey N. Weitzel, MD²⁴; Myra J. Wick, MD, PhD²⁵; Georgia L. Wiesner, MD, MS²⁶; Mary Dwyer, MS²⁷; and Susan Darlow, PhD²⁸

Abstract

The NCCN Clinical Practice Guidelines in Oncology for Genetic/Familial High-Risk Assessment: Breast and Ovarian provide recommendations for genetic testing and counseling for hereditary cancer syndromes and risk management recommendations for patients who are diagnosed with a syndrome. Guidelines focus on syndromes associated with an increased risk of breast and/or ovarian cancer. The NCCN Genetic/Familial High-Risk Assessment: Breast and Ovarian panel meets at least annually to review comments from reviewers within their institutions, examine relevant new data from publications and abstracts, and reevaluate and update their recommendations. The NCCN Guidelines Insights summarize the panel's discussion and most recent recommendations regarding risk management for carriers of moderately penetrant genetic mutations associated with breast and/or ovarian cancer.

J Natl Compr Canc Netw 2017;15(1):9-20

From ¹Fox Chase Cancer Center; ²The Ohio State University Comprehensive Cancer Center—James Cancer Hospital and Solow Research Institute; ³St. Jude Children's Research Hospital/The University of Tennessee Health Science Center; ⁴Huntsman Cancer Institute at the University of Utah; ⁵University of Alabama at Birmingham Comprehensive Cancer Center; ⁶FORCE: Facing Our Risk of Cancer Empowered; Dana-Farber/Boston and Women's Cancer Center; ⁷Duke Cancer Institute; ⁸Robert H. Lurie Comprehensive Cancer Center of Northwestern University; ⁹University of Colorado Cancer Center; ¹⁰Stanford Cancer Institute; ¹¹The University of Texas MD Anderson Cancer Center; ¹²UC San Diego Moores Cancer Center; ¹³University of Michigan Comprehensive Cancer Center; ¹⁴Memorial Sloan Kettering Cancer Center; ¹⁵Moffitt Cancer Center; ¹⁶Fred & Pamela Buffett Cancer Center; ¹⁷Massachusetts General Hospital Cancer Center; ¹⁸University of Washington Medical Center/Seattle Cancer Care Alliance; ¹⁹Case Comprehensive Cancer Center/University Hospitals Seidman Cancer Center and Cleveland Clinic Taussig Cancer Institute; ²⁰Roswell Park Cancer Institute; ²¹City of Hope Comprehensive Cancer Center; ²²Mayo Clinic Cancer Center; ²³Vanderbilt-Ingram Cancer Center; and ²⁴National Comprehensive Cancer Network.

Please Note

The NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) are a statement of consensus of the authors regarding their views of currently accepted approaches to treatment. The NCCN Guidelines® Insights highlight important changes to the NCCN Guidelines® recommendations from previous versions. Colored markings in the algorithm show changes and the discussion aims to further the understanding of these changes by summarizing salient portions of the NCCN Guideline Panel discussion, including the literature reviewed. These NCCN Guidelines Insights do not represent the full NCCN Guidelines; further, the National Comprehensive Cancer Network® (NCCN®) makes no representation or warranties of any kind regarding the content, use, or application of the NCCN Guidelines and NCCN Guidelines Insights and disclaims any responsibility for their applications or use in any way.

The full and most current version of these NCCN Guidelines are available at NCCN.org.

© National Comprehensive Cancer Network, Inc. 2017. All rights reserved. The NCCN Guidelines and the illustrations herein may not be reproduced in any form without the express written permission of NCCN.

*Provided content development and/or authorship assistance.

NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®)


Genetic/Familial High-Risk Assessment: Breast and Ovarian

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9

NCCN Guidelines® Insights Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

Featured Updates to the NCCN Guidelines

Mary B. Daly, MD, PhD^{1,2}; Robert Pilarski, MS, CGC³; Michael Berry, MD⁴; Sandra S. Bays, MD⁵; Meagan Farmer, MS, CGC⁶; Susan Friedman, DVM⁷; Judy E. Garber, MD, MPH⁸; Noah D. Kauff, MD⁹; Seema Khan, MD¹⁰; Catherine Klein, MD¹¹; Wendy Kohlmann, MS, CGC¹²; Allison Kurian, MD, MS¹³; Jennifer K. Litton, MD¹⁴; Lisa Madlensky, PhD, CGC¹⁵; Sofia D. Merajver, MD, PhD¹⁶; Kenneth Offit, MD¹⁷; Toyo Pal, MD¹⁸; Gwen Reiser, MS, CGC¹⁹; Kristen Mahoney Shannon, MS, CGC²⁰; Elizabeth Swisher, MD²¹; Shweta Vinayak, MD²²; Nicoleta C. Voian, MD, MPH²³; Jeffrey N. Weitzel, MD²⁴; Myra J. Wick, MD, PhD²⁵; Georgia L. Wiesner, MD, MS²⁶; Mary Dwyer, MS²⁷; and Susan Darlow, PhD²⁸

Abstract

The NCCN Clinical Practice Guidelines in Oncology for Genetic/Familial High-Risk Assessment: Breast and Ovarian provide recommendations for genetic testing and counseling for hereditary cancer syndromes and risk management recommendations for patients who are diagnosed with a syndrome. Guidelines focus on syndromes associated with an increased risk of breast and/or ovarian cancer. The NCCN Genetic/Familial High-Risk Assessment: Breast and Ovarian panel meets at least annually to review comments from reviewers within their institutions, examine relevant new data from publications and abstracts, and reevaluate and update their recommendations. The NCCN Guidelines Insights summarize the panel's discussion and most recent recommendations regarding risk management for carriers of moderately penetrant genetic mutations associated with breast and/or ovarian cancer.

J Natl Compr Canc Netw 2017;15(1):9-20

From ¹Fox Chase Cancer Center; ²The Ohio State University Comprehensive Cancer Center—James Cancer Hospital and Solow Research Institute; ³St. Jude Children's Research Hospital/The University of Tennessee Health Science Center; ⁴Huntsman Cancer Institute at the University of Utah; ⁵University of Alabama at Birmingham Comprehensive Cancer Center; ⁶FORCE: Facing Our Risk of Cancer Empowered; Dana-Farber/Boston and Women's Cancer Center; ⁷Duke Cancer Institute; ⁸Robert H. Lurie Comprehensive Cancer Center of Northwestern University; ⁹University of Colorado Cancer Center; ¹⁰Stanford Cancer Institute; ¹¹The University of Texas MD Anderson Cancer Center; ¹²UC San Diego Moores Cancer Center; ¹³University of Michigan Comprehensive Cancer Center; ¹⁴Memorial Sloan Kettering Cancer Center; ¹⁵Moffitt Cancer Center; ¹⁶Fred & Pamela Buffett Cancer Center; ¹⁷Massachusetts General Hospital Cancer Center; ¹⁸University of Washington Medical Center/Seattle Cancer Care Alliance; ¹⁹Case Comprehensive Cancer Center/University Hospitals Seidman Cancer Center and Cleveland Clinic Taussig Cancer Institute; ²⁰Roswell Park Cancer Institute; ²¹City of Hope Comprehensive Cancer Center; ²²Mayo Clinic Cancer Center; ²³Vanderbilt-Ingram Cancer Center; and ²⁴National Comprehensive Cancer Network.

*Provided content development and/or authorship assistance.

Please Note

The NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines[®]) are a statement of consensus of the authors regarding their views of currently accepted approaches to treatment. The NCCN Guidelines[®] Insights highlight important changes to the NCCN Guidelines[®] recommendations from previous versions. Colored markings in the algorithm show changes and the discussion aims to further the understanding of these changes by summarizing salient portions of the NCCN Guideline Panel discussion, including the literature reviewed. These NCCN Guidelines Insights do not represent the full NCCN Guidelines; further, the National Comprehensive Cancer Network[®] (NCCN[®]) makes no representation or warranties of any kind regarding the content, use, or application of the NCCN Guidelines and NCCN Guidelines Insights and disclaims any responsibility for their applications or use in any way.

The full and most current version of these NCCN Guidelines are available at NCCN.org.

© National Comprehensive Cancer Network, Inc. 2017. All rights reserved. The NCCN Guidelines and the illustrations herein may not be reproduced in any form without the express written permission of NCCN.

NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines[®])

Genetic/Familial High-Risk Assessment: Breast and Ovarian

Version 2.2016

NCCN.org

Continue

BRCA+16GENES: International Guidelines

Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

BREAST AND OVARIAN MANAGEMENT BASED ON GENETIC TEST RESULTS^{a,b}

The inclusion of a gene on this table below does not imply the endorsement either for or against multi-gene testing for moderate-penetrance genes.

Gene	Breast Cancer Risk and Management	Ovarian Cancer Risk and Management	Other Cancer Risks and Management
<i>ATM</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast starting at age 40 y^c RRM: Consider based on family history 	No increased risk of OC	Unknown or insufficient evidence for pancreas or prostate cancer
	Comments: Insufficient evidence to recommend against radiation therapy. The 7271T>G missense mutation may act in a dominant-negative fashion, resulting in a lifetime breast cancer risk as high as 60% by age 80 (which is higher than truncating mutations, where risks are in the range of 30-40%). Counsel for risk of autosomal recessive condition in offspring.		
<i>BRCA1</i>	Increased risk of BC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Increased risk of OC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Prostate cancer <ul style="list-style-type: none"> See BRCA Mutation-Positive Management
<i>BRCA2</i>	Increased risk of BC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Increased risk of OC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Pancreas, Prostate, Melanoma <ul style="list-style-type: none"> See BRCA Mutation-Positive Management
<i>BRIP1</i>	No increased risk of BC	Increased risk of OC <ul style="list-style-type: none"> Consider RRSO at 45–50 y 	N/A
	Comments: Counsel for risk of autosomal recessive condition in offspring. Based on estimates from available studies, the lifetime risk of ovarian cancer in carriers of mutations in <i>BRIP1</i> appears to be sufficient to justify consideration of risk-reducing salpingo-oophorectomy. The current evidence is insufficient to make a firm recommendation as to the optimal age for this procedure. Based on the current, limited evidence base, a discussion about surgery should be held around age 45–50 y or earlier based on a specific family history of an earlier onset ovarian cancer.		
<i>CDH1</i>	Increased risk of lobular BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast starting at age 30 y^c RRM: Consider based on family history 	No increased risk of OC	Diffuse gastric cancer <ul style="list-style-type: none"> See NCCN Guidelines for Gastric Cancer

BC: Breast cancer
OC: Ovarian cancer
RRM: Risk-reducing mastectomy
RRSO: Risk-reducing salpingo-oophorectomy

^aTung N, Domchek SM, Stadler Z, Nathanson KL, Couch F, Garber JE, Offit K, Robson ME. Counselling framework for moderate-penetrance cancer-susceptibility mutations. *Nat Rev Clin Oncol* 2016;13:581-588.

^bThe following genes and others are found on some of the panels but there is insufficient evidence to make any recommendations for breast MRI, RRSO, or RRM: BARD1, FANCC, MRE11A, MUTYH heterozygotes, REQL, RAD50, RET1, SLX4, SMARCA4, or XRCC2.

^cMay be modified based on family history or specific gene mutation.

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Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

BREAST AND OVARIAN MANAGEMENT BASED ON GENETIC TEST RESULTS^{a,b}

The inclusion of a gene on this table below does not imply the endorsement either for or against multi-gene testing for moderate-penetrance genes.

Gene	Breast Cancer Risk and Management	Ovarian Cancer Risk and Management	Other Cancer Risks and Management
<i>ATM</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast starting at age 40 y^c RRM: Consider based on family history 	No increased risk of OC	Unknown or insufficient evidence for pancreas or prostate cancer
	Comments: Insufficient evidence to recommend against radiation therapy. The 7271T>G missense mutation may act in a dominant-negative fashion, resulting in a lifetime breast cancer risk as high as 60% by age 80 (which is higher than truncating mutations, where risks are in the range of 30-40%). Counsel for risk of autosomal recessive condition in offspring.		
<i>BRCA1</i>	Increased risk of BC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Increased risk of OC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Prostate cancer <ul style="list-style-type: none"> See BRCA Mutation-Positive Management
<i>BRCA2</i>	Increased risk of BC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Increased risk of OC <ul style="list-style-type: none"> See BRCA Mutation-Positive Management 	Pancreas, Prostate, Melanoma <ul style="list-style-type: none"> See BRCA Mutation-Positive Management
<i>BRIP1</i>	No increased risk of BC	Increased risk of OC <ul style="list-style-type: none"> Consider RRSO at 45–50 y 	N/A
	Comments: Counsel for risk of autosomal recessive condition in offspring. Based on estimates from available studies, the lifetime risk of ovarian cancer in carriers of mutations in <i>BRIP1</i> appears to be sufficient to justify consideration of risk-reducing salpingo-oophorectomy. The current evidence is insufficient to make a firm recommendation as to the optimal age for this procedure. Based on the current, limited evidence base, a discussion about surgery should be held around age 45–50 y or earlier based on a specific family history of an earlier onset ovarian cancer.		
<i>CDH1</i>	Increased risk of lobular BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast starting at age 30 y^c RRM: Consider based on family history 	No increased risk of OC	Diffuse gastric cancer <ul style="list-style-type: none"> See NCCN Guidelines for Gastric Cancer

BC: Breast cancer
OC: Ovarian cancer
RRM: Risk-reducing mastectomy
RRSO: Risk-reducing salpingo-oophorectomy

^aTung N, Domchek SM, Stadler Z, Nathanson KL, Couch F, Garber JE, Offit K, Robson ME. Counselling framework for moderate-penetrance cancer-susceptibility mutations. *Nat Rev Clin Oncol* 2016;13:581-588.

^bThe following genes and others are found on some of the panels but there is insufficient evidence to make any recommendations for breast MRI, RRSO, or RRM: BARD1, FANCC, MRE11A, MUTYH heterozygotes, REQL, RAD50, RET1, SLX4, SMARCA4, or XRCC2.

^cMay be modified based on family history or specific gene mutation.

BRCA+16GENES: International Guidelines

Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

BREAST AND OVARIAN MANAGEMENT BASED ON GENETIC TEST RESULTS^a

The inclusion of a gene on this table below does not imply the endorsement either for or against multi-gene testing for moderate-penetrance genes.

Gene	Breast Cancer Risk and Management	Ovarian Cancer Risk and Management	Other Cancer Risks and Management
<i>CHEK2</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast age 40 y^c RRM: Evidence insufficient, manage based on family history. 	No increased risk of OC	Colon <ul style="list-style-type: none"> See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal
Comments: Risk data are based only on frameshift mutations. The risks for most missense mutations are unclear.			
<i>MSH2, MLH1, MSH6, PMS2, EPCAM</i>	Unknown or insufficient evidence for BC risk^d <ul style="list-style-type: none"> Manage based on family history 	Increased risk of OC <ul style="list-style-type: none"> See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal 	See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal
<i>NBN</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast age 40 y^c RRM: Evidence insufficient, manage based on family history 	Unknown or insufficient evidence for OC risk	Unknown or insufficient evidence
Comments: Management recommendations are based on data derived from the 657del5 Slavic truncating mutation. Although risks for other mutations have not been established it is prudent to manage patients with other truncating mutations similarly to those with 675del5. Counsel for risk of autosomal recessive condition in children.			
<i>NF1</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram starting at age 30 y and consider breast MRI with contrast from ages 30–50 y RRM: Evidence insufficient, manage based on family history. 	No increased risk of OC	<ul style="list-style-type: none"> Malignant peripheral nerve sheath tumors, GIST, others Recommend referral to NF specialist for evaluation and management.
Comments: At this time, there are no data to suggest an increased breast cancer risk after age 50 y.			

^aTung N, Domchek SM, Stadler Z, Nathanson KL, Couch F, Garber JE, Offit K, Robson ME. Counselling framework for moderate-penetrance cancer-susceptibility mutations. *Nat Rev Clin Oncol* 2016;13:581-588.

^cMay be modified based on family history or specific gene mutation.

^dThere have been suggestions that there is an increased risk for breast cancer in LS patients; however, there is not enough evidence to support increased screening above average-risk breast cancer screening recommendations.

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Genetic/Familial High-Risk Assessment: Breast and Ovarian, Version 2.2017

BREAST AND OVARIAN MANAGEMENT BASED ON GENETIC TEST RESULTS^a

The inclusion of a gene on this table below does not imply the endorsement either for or against multi-gene testing for moderate-penetrance genes.

Gene	Breast Cancer Risk and Management	Ovarian Cancer Risk and Management	Other Cancer Risks and Management
<i>CHEK2</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast age 40 y^c RRM: Evidence insufficient, manage based on family history. 	No increased risk of OC	Colon <ul style="list-style-type: none"> See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal
Comments: Risk data are based only on frameshift mutations. The risks for most missense mutations are unclear.			
<i>MSH2, MLH1, MSH6, PMS2, EPCAM</i>	Unknown or insufficient evidence for BC risk^d <ul style="list-style-type: none"> Manage based on family history 	Increased risk of OC <ul style="list-style-type: none"> See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal 	See NCCN Guidelines for Genetic/Familial High-Risk Assessment: Colorectal
<i>NBN</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram and consider breast MRI with contrast age 40 y^c RRM: Evidence insufficient, manage based on family history 	Unknown or insufficient evidence for OC risk	Unknown or insufficient evidence
Comments: Management recommendations are based on data derived from the 657del5 Slavic truncating mutation. Although risks for other mutations have not been established it is prudent to manage patients with other truncating mutations similarly to those with 675del5. Counsel for risk of autosomal recessive condition in children.			
<i>NF1</i>	Increased risk of BC <ul style="list-style-type: none"> Screening: Annual mammogram starting at age 30 y and consider breast MRI with contrast from ages 30–50 y RRM: Evidence insufficient, manage based on family history. 	No increased risk of OC	<ul style="list-style-type: none"> Malignant peripheral nerve sheath tumors, GIST, others Recommend referral to NF specialist for evaluation and management.
Comments: At this time, there are no data to suggest an increased breast cancer risk after age 50 y.			

^aTung N, Domchek SM, Stadler Z, Nathanson KL, Couch F, Garber JE, Offit K, Robson ME. Counselling framework for moderate-penetrance cancer-susceptibility mutations. Nat Rev Clin Oncol 2016;13:581-588.

^cMay be modified based on family history or specific gene mutation.

^dThere have been suggestions that there is an increased risk for breast cancer in LS patients; however, there is not enough evidence to support increased screening above average-risk breast cancer screening recommendations.

BC: Breast cancer

OC: Ovarian cancer

RRM: Risk-reducing mastectomy

RRSO: Risk-reducing salpingo-oophorectomy

Continued

BRCA+16GENES: BRCA mutation patient management

BRCA mutations

1. Breast awareness starting at age **18 y**.
2. Clinical breast exam, every 6–12 mo, starting at age **25 y**.
3. Breast screening
 - Age **25–29 y**, annual breast MRI screening (preferred) or mammogram if MRI is unavailable or individualized based on family history if a breast cancer diagnosis before age 30 is present.
 - Age **30–75 y**, annual mammogram and breast MRI screening.
 - Age **>75 y**, management should be considered on an individual basis.
4. For women with a BRCA mutation who are treated for breast cancer, screening of remaining breast tissue with annual mammography and breast MRI should continue.
5. Discuss option of risk-reducing mastectomy. Counselling may include a discussion regarding degree of protection, reconstruction options, and risks.
6. Recommend risk-reducing salpingo-oophorectomy (RRSO), typically between 35 and 40 y, and upon completion of child bearing. Because ovarian cancer onset in patients with *BRCA2* mutations is an average of 8–10 years later than in patients with *BRCA1* mutations, it is reasonable to delay RRSO until age 40–45 y in patients with *BRCA2* mutations who have already maximized their breast cancer prevention (i.e., undergone bilateral mastectomy).

CDH1 mutations

Germline mutations in *CDH1* have reported a **cumulative lifetime risk for breast cancer of 39% to 52%** NCCN recommends **screening with annual mammogram** (or consideration of breast MRI) **beginning at age 30 years**. Screening may be considered earlier in patients with a family history of early-onset breast cancer. The option of **risk-reducing mastectomy should be discussed for these carriers**.

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
BRCA+16GENES: Analysis technology


Next-generation sequencing (NGS) with paired-end reads

- The sequencing process consists of several steps. First, samples must be prepared for sequencing. To do this, DNA is extracted, purified and fragmented, and adapters are added.
- The sample is placed on a glass slide called a flow cell, which acts as a working surface. The surface of the flow cell features additional probes with adapters added to the DNA so that the sample may be fixed to the surface. Next, sequencing by synthesis begins. This process generates large amounts of information, which will be analysed using a sophisticated software system.

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
Next-generation sequencing (NGS) με paired-end reads

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BRCA+16GENES: Analysis technology


MLPA (*Multiplex Ligation-dependent Probe Amplification*)

- Performing the MLPA technique requires carrying out an initial DNA denaturing process and then a hybridisation process with specific probes. Once the probes have bound to the study region, the products are simultaneously amplified by PCR.
- The products are analysed using electrophoresis, which separates the fragments by size in order to detect potential deletions and duplications.

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BRCA+16GENES: Indications from various medical associations

Patients suffering from cancer					
INDICATION	SEOM	SGO	ASCO	ACMG	ESMO
Women with high-grade epithelial carcinoma regardless of age		●	●	●	
Breast cancer at or before 50 years of age				●	
Breast cancer at or before 45 years of age		●	●		
Breast cancer at or before 40 years of age	●				
Bilateral breast cancer or diagnosis of two primary breast cancers, the first being before 50 years of age	●	●	●	●	●
Triple-negative breast cancer before 60 years of age		●		●	●
Breast cancer with first-, second- or third-degree relative with breast cancer at or before 50 years of age or ovarian cancer at any age	●	●	●	●	
Breast cancer with two or more relatives with pancreatic cancer or prostate cancer (Gleason > 7)		●		●	
Breast cancer in males	●	●	●		●

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ASCO: American Society of Clinical Oncology
ACMG: American College of Medical Genetics and Genomics
ESMO: European Society for Medical Oncology

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ASCO: American Society of Clinical Oncology
 ACMG: American College of Medical Genetics and Genomics
 ESMO: European Society for Medical Oncology

BRCA+16GENES: Indications from various medical associations

Relatives of patients suffering from cancer				
INDICATION	SEOM	SGO	ASCO	ACMG
Person with a male relative with breast cancer	●	●	●	
Person with a relative with breast cancer and another relative with diffuse gastric cancer, at least one with a diagnosis before 50 years of age				●
Person with a relative who is a carrier of a known mutation in <i>BRCA1/BRCA2</i>	●	●	●	●
Person with a first-degree relative with breast cancer before 45 years of age or with ovarian cancer at any age	●	●	●	●
Person with a first-degree relative or more than one first-, second- or third-degree relative with breast cancer and descendent from Ashkenazi Jews		●		
Person with a first-degree relative or more than one first-, second- or third-degree relative with triple-negative breast cancer before 60 years of age		●		

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ASCO: American Society of Clinical Oncology
ACMG: American College of Medical Genetics and Genomics

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BRCA+16GENES: Indications

- Women with a **family history of breast cancer** (female or male) and/or **ovarian cancer**.
- Women with a **relative affected with family cancer syndrome**.
- **Patients with these type of tumours** in order to determine their potential hereditary nature.
- **Women ≥ 30 years of age with no prior family history**, to determine the genetic risk of breast and ovarian hereditary cancer and evaluate the different preventive options.¹

First-degree relatives of a carrier of the mutation have a higher risk of being a carrier and developing the disease.

1. King MC, Levy-Lahad E, Lahad A. Population-Based Screening for BRCA1 and BRCA2: 2014 Lasker Award. JAMA. 2014;312(11):1091-2.

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
BRCA+16GENES: Important considerations

The presence of mutations in the genes of the panel associated to breast, ovarian and endometrial cancer may also imply a **higher risk for other type of cancers or hereditary cancer syndromes (e.g. Lynch, Cowden and Li-Fraumeni syndromes)**. This information, if applies, may be included in the results report.

SYNLAB **BRCA+16 GENES**

Patient report
Hereditary predisposition to breast, ovarian and endometrial cancer

PATIENT INFORMATION	ORDERING PHYSICIAN/CLINIC INFORMATION
First name:	Doctor:
Last names:	Medical speciality:
Gender:	Center/Hospital:
Date of birth (Day/Month/Year):	

SAMPLE INFORMATION	
Draw date:	Sample ID:
Receipt date:	Sample origin ID:
	Report date:

REASON FOR THE STUDY

RESULT

NO MUTATION DETECTED

INTERPRETATION


Tip: here
No mutations classified as pathogenic, probably pathogenic or of uncertain significance have been identified in this study. A result without pathogenic or probably pathogenic mutation detected does not eliminate the risk of breast, ovary and endometrial cancer.
This report is to be interpreted in the clinical context and family history of the patient.

SYNLAB DIAGNOSTICOS GLOBALES SL.
C/ Valgrande 8, Edificio Tharwerrh 3, 28108 Alcobendas, Madrid, España.
www.synlab.com / 900 400 442

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SYNLAB  **BRCA+16 GENES**

Patient report
Hereditary predisposition to breast, ovarian and endometrial cancer

PATIENT INFORMATION	ORDERING PHYSICIAN/CLINIC INFORMATION
First name:	Doctor:
Last names:	Medical speciality:
Gender:	Center/Hospital:
Date of birth (Day/Month/Year):	

SAMPLE INFORMATION
Draw date:
Receipt date:
Sample ID:
Sample origin ID:
Report date:

REASON FOR THE STUDY

RESULT

NO MUTATION DETECTED

INTERPRETATION

Tipos here
No mutations classified as pathogenic, probably pathogenic or of uncertain significance have been identified in this study. A result without pathogenic or probably pathogenic mutation detected does not eliminate the risk of breast, ovary and endometrial cancer.
This report is to be interpreted in the clinical context and family history of the patient.

SYNLAB DIAGNOSTICOS GLOBALES SL.
C/ Valgrande 8, Edificio Tharwenth 3, 28108 Alcobendas, Madrid, España.
www.synlab.com / 900 400 442

Page 1 of 3

BRCA+16GENES: Advantages

✓ COMPREHENSIVE

Includes the genes with robust scientific evidence related with these types of cancer, for which **a specific patient management is described on the NCCN guidelines**, not solely *BRCA1* and *BRCA2*.

NCCN: National Comprehensive Cancer Network.

✓ RELIABLE

- NGS sequencing with paired-end reads of the genes included on the panel.
- Duplications and deletions in *BRCA1*, *BRCA2* and *EPCAM* genes.
- Positive results confirmation with Sanger sequencing.

✓ EASY AND SIMPLE


Saliva or blood sample available.

✓ DATA BASES

Classification and study of variants with the most complete databases.

✓ COMPREHENSIVE AND SIMPLE TEST REPORT

Report designed to facilitate the result interpretation.

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
BRCA+16GENES: Advantages

✓ **FAST AND AFORDABLE**

The high degree of automation allows to obtain results in 10 working days.

✓ **PERFORMED ENTIRELY AT SYNLAB'S LABORATORIES**

✓ **SYNLAB QUALITY AND EXPERTISE**

Developed by the genetic experts of SYNLAB group, Europe's number one medical diagnostics provider.

✓ **GENETIC COUNSELLING**

SYNLAB put at your disposal without additional cost:

With access to our genetic counselling platform where you will receive genetic advice from our experts via videoconference. You can also contact us sending an e-mail to:
genetic.counselling@synlab.com

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BRCA+16GENES: Sampling kit

BRCA+16 GENES kit

The sample collection kit includes:

- Device to collect saliva
- Test requisition form and informed consent
- Instructions
- Material for sending the sample

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BRCA+16GENES: Test Requisition Form and Informed Consent

BRCA+16 GENES

Test Requisition Form

SYNLAB BARCODE LABEL

PATIENT INFORMATION							
Name:		Surname(s):					
Gender: <input type="checkbox"/> Male <input type="checkbox"/> Female	National ID #:	Date of birth: / / (day/month/year)					
Telephone:	Email:	Address:					
SAMPLE INFORMATION							
Type of sample: <input type="checkbox"/> Blood <input type="checkbox"/> Saliva	Sample ID:	Draw date: / / (day/month/year)					
ORDERING PHYSICIAN / CENTRE INFORMATION							
Name and surname(s):		Medical license #:	Client code:				
Email:		Telephone:	Specialty:				
Centre/hospital name:		Address:					
PATIENT'S CLINICAL HISTORY							
Ancestry (tick all that apply):							
<input type="checkbox"/> Western/Northern European	<input type="checkbox"/> Central/Eastern European	<input type="checkbox"/> Southern European	<input type="checkbox"/> African				
<input type="checkbox"/> Latin American/Caribbean	<input type="checkbox"/> Native American	<input type="checkbox"/> Near Eastern/Middle Eastern	<input type="checkbox"/> Asian				
<input type="checkbox"/> Ashkenazi	<input type="checkbox"/> Other: _____						
Patient's personal history of cancer (tick all that apply):							
<input type="checkbox"/> No personal history of cancer							
<input type="checkbox"/> Breast cancer - age at diagnosis: _____ years	<input type="checkbox"/> Bilateral	<input type="checkbox"/> Premenopausal	<input type="checkbox"/> Triple-negative				
<input type="checkbox"/> Ovarian cancer - age at diagnosis: _____ years							
<input type="checkbox"/> Other type: _____	- age at diagnosis: _____ years						
Family history of cancer:							
<input type="checkbox"/> No known family history <input type="checkbox"/> Known family history - fill in the table below:							
If breast cancer							
Type of cancer	Family relationship	Maternal	Paternal	Age at diagnosis	Bilateral	Premenopausal	Triple-negative
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Other significant information:							
<ul style="list-style-type: none"> • Smoker: <input type="checkbox"/> Yes <input type="checkbox"/> No • Has received a bone marrow transplant: <input type="checkbox"/> Yes <input type="checkbox"/> No • In women who have given birth, breast-feeding for at least 1 year: <input type="checkbox"/> Yes <input type="checkbox"/> No 							
ORDERING PHYSICIAN'S SIGNATURE							
By signing this form I certify that, prior to performing the BRCA+16GENES test, I have informed the patient of the risks and implications that performing this test represents. I certify that all the patient's questions have been resolved and that I have received the patient's explicit consent to perform the test.							
Physician's signature:		Draw: / / (day/month/year)					
PATIENT'S INFORMED CONSENT							
By signing this form I confirm that I have read and agreed to the information contained on both sides of the form, or that it has been read to me, and that I have understood its content. I have received genetic counselling from my physician (or another person indicated by my physician) regarding the purpose of the test and its potential risks and limitations. I have been given the opportunity to ask any questions I had. I have received answers to all my questions and I have been given enough time to reflect on the information and my decision to undergo this test. I consent to undergo this test and discuss the results and appropriate medical management with my specialist. I agree to my biological sample being used solely and exclusively for the test specified in this requisition form, and for no other type of test under any circumstances. I understand that my physician has determined that this test is appropriate for me and I authorise SYNLAB DIAGNOSTICOS GLOBALES SA with corporate tax ID no. A-59845875, and registered address at C/ Verge de Guadalupe 18, 08550 Espoluges de Llobregat, Spain, [INSERT LOCAL SYNLAB COUNTRY ENTITY] with registered office at [INSERT LOCAL SYNLAB COUNTRY ENTITY ADDRESS] and its affiliates (collectively referred to as "SYNLAB") to perform the BRCA+16GENES test, as well as the results being sent to my specialist. By signing below I agree to the foregoing and the terms of the Patient Informed Consent.							
Patient's or legal representative's signature:		Draw: / / (day/month/year)					
BILLING DETAILS							
<input type="checkbox"/> Patient <input type="checkbox"/> Insurance company: _____ <input type="checkbox"/> Ordering physician <input type="checkbox"/> Other: _____							

BRCA+16 GENES

Patient Informed Consent

(Copy for the laboratory)

The BRCA+16GENES test consists of extracting and quantifying the DNA in the sample received and DNA sequencing by Next-Generation Sequencing (NGS). This technique is capable of detecting point mutations and small insertions/deletions throughout the coding sequence and the flanking intronic region of the ATM, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, EPCAM, MLH1, MSH2, MSH6, NBN, PALB2, PMS2, PTEN, RAD51C, RAD51D, STK11 and TP53 genes. The technique is supplemented with Multiplex Ligation-dependent Probe Amplification (MLPA) of the BRCA1, BRCA2 and EPCAM genes to detect large deletions and duplications. The pathogenic and probably pathogenic variants detected using next-generation sequencing are confirmed through Sanger sequencing.

The genes included in the BRCA+16GENES test are involved in cell cycle control and DNA repair during cell division. Mutations in these genes lead to a loss of cell control and capacity for DNA repair, which may imply a greater risk of developing cancer than that of the general population. Abnormalities in the genes included in the BRCA+16GENES test represent an increased risk of suffering from hereditary breast, ovarian and endometrial cancer. Mutations can be inherited from both the mother and the father. The probability of transmitting a mutation to the offspring is 50%. The presence of a *de novo* (non-inherited) mutation cannot be excluded.

The BRCA+16GENES test has certain limitations in determining the risk of the patient and/or the patient's family members of having hereditary breast, ovarian and endometrial cancer. Identifying a pathological genetic abnormality represents an increased risk of having the associated disease, but does not necessarily imply its development. If a relevant genetic abnormality is identified, it will be useful to perform a genetic study of the immediate family members (parents, children, siblings, etc.). If the family members tested did not have the previously detected abnormality, this would mean that the risk of eventually developing the disease is not increased, that is to say, is equal to that of the general population. Should pathological variants not be detected, this DOES NOT ELIMINATE the possibility of the patient having cancer or other genetic diseases. The presence of mutations in the genes of the panel related to breast, ovarian and endometrial cancer may also imply a higher risk for other cancer types or hereditary cancer syndromes (e. g. Lynch, Cowden and Li-Fraumeni syndromes). This information, if applicable, will be included on the report.

In compliance with the provisions of the current legislation, the patient accepts and recognizes at all the effects that the service will be deemed to have been fully executed once the sample has been taken. Once such circumstance has occurred, the patient will lose its right to cancel the contracted service and SYNLAB will not be obliged to reimburse the amounts received for that service.

The sample will be stored as set out by the regulations that apply to clinical diagnostic laboratories. Once the result has been issued, if there is any surplus sample, it will be stored for 1 month after being analysed, and once this period has elapsed it will be destroyed. The result of the BRCA+16GENES test is confidential. The patient's results will only be given to the patient's physician or another professional involved in the patient's medical care, unless the communication of this information is required by mandatory law or by order of authorities to disclose the patient's data to authorities or third parties, such as regulatory authorities. The healthcare professional is responsible for explaining the specific use and limitations of this test to the patient. It is recommended that the results are reported to the patient by a specialist in a medical consultation. The result may occasionally be delayed or a second sample required. Performing the test does not include free direct genetic counselling for the patient; however, SYNLAB offers this service through its Genetic Counselling Unit (genetic.counselling@synlab.com).

The patient agrees that their biological samples and a copy of this executed Test Requisition Form and all personal data about them contained in this form are transferred to and processed by the laboratory SYNLAB DIAGNOSTICOS GLOBALES SA with registered office at C/ Verge de Guadalupe 18, 08550 Espoluges de Llobregat, Spain, and that the test results and the personal data may be processed and stored both by [INSERT LOCAL SYNLAB COUNTRY ENTITY] with registered office at [INSERT LOCAL SYNLAB COUNTRY ENTITY ADDRESS] and its affiliates (collectively referred to as "SYNLAB"), where the level of protection may not be the same as in the patient's country. The patient has a right to withdraw their consent, but in this case it is understood that the BRCA+16GENES test cannot be provided. The patient's personal data will only be used to perform the test, to communicate with them and for invoicing purposes. The patient further understands that their personal data will be stored for a duration of [] after the test has been performed; that they may exercise the rights of access, rectification and, as applicable, restriction, opposition or erasure by sending an email to [INSERT LOCAL SYNLAB COUNTRY ENTITY AND ADDRESS] at [INSERT LOCAL SYNLAB COUNTRY ENTITY ADDRESS] and that they have a right to lodge a complaint with the competent Supervisory Authority in their country. They also agree that the results of the test will be communicated by the testing laboratory to [INSERT LOCAL SYNLAB COUNTRY ENTITY] and to the physician mentioned on this form or its/their representatives.

Pursuant to the best practices and quality standards of clinical laboratories, the patient acknowledges that SYNLAB may use the leftover specimen and the patient's medical and genetic information, in an anonymized form (unless forbidden by applicable legislation) for research or quality assurance purposes. Such uses may result in the development of commercial products and services. The patient will not receive notice of any specific uses and will not receive any compensation for these uses. In any event, all such uses will be in compliance with applicable legislation.

Tick the box if you would not like your sample to be used for research purposes.

In compliance with Law 14/2007, of 3 July, on Biomedical Research, specifically its Articles 47 and 48, the prescribing physician must obtain informed consent to perform genetic tests. The patient's signature on this consent form is aimed at complying with this requirement.

Patient's name and surname(s): _____

Patient's or legal representative's signature: _____ Date: / / (day/month/year)

TRF INT ENG_V1 900 400 442 • atencion.cliente@synlab.com • www.synlab.com pag. 2 / 4

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Test Requisition Form

SYNLAB BARCODE LABEL

PATIENT INFORMATION							
Name:		Surname(s):					
Gender: <input type="checkbox"/> Male <input type="checkbox"/> Female	National ID #:	Date of birth: / / (day/month/year)					
Telephone:	Email:	Address:					
SAMPLE INFORMATION							
Type of sample: <input type="checkbox"/> Blood <input type="checkbox"/> Saliva	Sample ID:	Draw date: / / (day/month/year)					
ORDERING PHYSICIAN / CENTRE INFORMATION							
Name and surname(s):		Medical license #:	Client code:				
Email:		Telephone:	Specialty:				
Centre/hospital name:		Address:					
PATIENT'S CLINICAL HISTORY							
Ancestry (tick all that apply):							
<input type="checkbox"/> Western/Northern European	<input type="checkbox"/> Central/Eastern European	<input type="checkbox"/> Southern European	<input type="checkbox"/> African				
<input type="checkbox"/> Latin American/Caribbean	<input type="checkbox"/> Native American	<input type="checkbox"/> Near Eastern/Middle Eastern	<input type="checkbox"/> Asian				
<input type="checkbox"/> Ashkenazi	<input type="checkbox"/> Other: _____						
Patient's personal history of cancer (tick all that apply):							
<input type="checkbox"/> No personal history of cancer							
<input type="checkbox"/> Breast cancer - age at diagnosis: _____ years	<input type="checkbox"/> Bilateral	<input type="checkbox"/> Premenopausal	<input type="checkbox"/> Triple-negative				
<input type="checkbox"/> Ovarian cancer - age at diagnosis: _____ years							
<input type="checkbox"/> Other type: _____	- age at diagnosis: _____ years						
Family history of cancer:							
<input type="checkbox"/> No known family history <input type="checkbox"/> Known family history - fill in the table below:							
If breast cancer							
Type of cancer	Family relationship	Maternal	Paternal	Age at diagnosis	Bilateral	Premenopausal	Triple-negative
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Other significant information:							
<ul style="list-style-type: none"> • Smoker: <input type="checkbox"/> Yes <input type="checkbox"/> No • Has received a bone marrow transplant: <input type="checkbox"/> Yes <input type="checkbox"/> No • In women who have given birth, breast-feeding for at least 1 year: <input type="checkbox"/> Yes <input type="checkbox"/> No 							
ORDERING PHYSICIAN'S SIGNATURE							
By signing this form I certify that, prior to performing the BRCA+16GENES test, I have informed the patient of the risks and implications that performing this test represents. I certify that all the patient's questions have been resolved and that I have received the patient's explicit consent to perform the test.							
Physician's signature:		Draw: / / (day/month/year)					
PATIENT'S INFORMED CONSENT							
By signing this form I confirm that I have read and agreed to the information contained on both sides of the form, or that it has been read to me, and that I have understood its content. I have received genetic counselling from my physician (or another person indicated by my physician) regarding the purpose of the test and its potential risks and limitations. I have been given the opportunity to ask any questions I had. I have received answers to all my questions and I have been given enough time to reflect on the information and my decision to undergo this test. I consent to undergo this test and discuss the results and appropriate medical management with my specialist. I agree to my biological sample being used solely and exclusively for the test specified in this requisition form, and for no other type of test under any circumstances. I understand that my physician has determined that this test is appropriate for me and I authorise SYNLAB DIAGNOSTICOS GLOBALES SA with corporate tax ID no. A-59945875, and registered address at C/ Verge de Guadalupe 18, 08550 Espoluges de Llobregat, Spain, [INSERT LOCAL SYNLAB COUNTRY ENTITY] with registered office at [INSERT LOCAL SYNLAB COUNTRY ENTITY ADDRESS] and its affiliates (collectively referred to as "SYNLAB") to perform the BRCA+16GENES test, as well as the results being sent to my specialist. By signing below I agree to the foregoing and the terms of the Patient Informed Consent.							
Patient's or legal representative's signature:		Draw: / / (day/month/year)					
BILLING DETAILS							
<input type="checkbox"/> Patient <input type="checkbox"/> Insurance company: _____ <input type="checkbox"/> Ordering physician <input type="checkbox"/> Other: _____							

BRCA+16 GENES

Patient Informed Consent

(Copy for the laboratory)

The BRCA+16GENES test consists of extracting and quantifying the DNA in the sample received and DNA sequencing by Next-Generation Sequencing (NGS). This technique is capable of detecting point mutations and small insertions/deletions throughout the coding sequence and the flanking intronic region of the ATM, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, EPCAM, MLH1, MSH2, MSH6, NBN, PALB2, PMS2, PTEN, RAD51C, RAD51D, STK11 and TP53 genes. The technique is supplemented with Multiplex Ligation-dependent Probe Amplification (MLPA) of the BRCA1, BRCA2 and EPCAM genes to detect large deletions and duplications. The pathogenic and probably pathogenic variants detected using next-generation sequencing are confirmed through Sanger sequencing.

The genes included in the BRCA+16GENES test are involved in cell cycle control and DNA repair during cell division. Mutations in these genes lead to a loss of cell control and capacity for DNA repair, which may imply a greater risk of developing cancer than that of the general population. Abnormalities in the genes included in the BRCA+16GENES test represent an increased risk of suffering from hereditary breast, ovarian and endometrial cancer. Mutations can be inherited from both the mother and the father. The probability of transmitting a mutation to the offspring is 50%. The presence of a *de novo* (non-inherited) mutation cannot be excluded.

The BRCA+16GENES test has certain limitations in determining the risk of the patient and/or the patient's family members of having hereditary breast, ovarian and endometrial cancer. Identifying a pathological genetic abnormality represents an increased risk of having the associated disease, but does not necessarily imply its development. If a relevant genetic abnormality is identified, it will be useful to perform a genetic study of the immediate family members (parents, children, siblings, etc.). If the family members tested did not have the previously detected abnormality, this would mean that the risk of eventually developing the disease is not increased, that is to say, is equal to that of the general population. Should pathological variants not be detected, this DOES NOT ELIMINATE the possibility of the patient having cancer or other genetic diseases. The presence of mutations in the genes of the panel related to breast, ovarian and endometrial cancer may also imply a higher risk for other cancer types or hereditary cancer syndromes (e. g. Lynch, Cowden and Li-Fraumeni syndromes). This information, if applicable, will be included on the report.

In compliance with the provisions of the current legislation, the patient accepts and recognizes at all the effects that the service will be deemed to have been fully executed once the sample has been taken. Once such circumstance has occurred, the patient will lose its right to cancel the contracted service and SYNLAB will not be obliged to reimburse the amounts received for that service.

The sample will be stored as set out by the regulations that apply to clinical diagnostic laboratories. Once the result has been issued, if there is any surplus sample, it will be stored for 1 month after being analysed, and once this period has elapsed it will be destroyed. The result of the BRCA+16GENES test is confidential. The patient's results will only be given to the patient's physician or another professional involved in the patient's medical care, unless the communication of this information is required by mandatory law or by order of authorities to disclose the patient's data to authorities or third parties, such as regulatory authorities. The healthcare professional is responsible for explaining the specific use and limitations of this test to the patient. It is recommended that the results are reported to the patient by a specialist in a medical consultation. The result may occasionally be delayed or a second sample required. Performing the test does not include free direct genetic counselling for the patient; however, SYNLAB offers this service through its Genetic Counselling Unit (genetic.counselling@synlab.com).

The patient agrees that their biological samples and a copy of this executed Test Requisition Form and all personal data about them contained in this form are transferred to and processed by the laboratory SYNLAB DIAGNOSTICOS GLOBALES SA with registered office at C/ Verge de Guadalupe 18, 08550 Espoluges de Llobregat, Spain, and that the test results and the personal data may be processed and stored both by [INSERT LOCAL SYNLAB COUNTRY ENTITY] with registered office at [INSERT LOCAL SYNLAB COUNTRY ENTITY ADDRESS] and its affiliates (collectively referred to as "SYNLAB"), where the level of protection may not be the same as in the patient's country. The patient has a right to withdraw their consent, but in this case it is understood that the BRCA+16GENES test cannot be provided. The patient's personal data will only be used to perform the test, to communicate with them and for invoicing purposes. The patient further understands that their personal data will be stored for a duration of [] after the test has been performed; that they may exercise the rights of access, rectification and, as applicable, restriction, opposition or erasure by sending an email to [INSERT LOCAL SYNLAB COUNTRY ENTITY AND ADDRESS] at [INSERT LOCAL SYNLAB COUNTRY ENTITY EMAIL], and that they have a right to lodge a complaint with the competent Supervisory Authority in their country. They also agree that the results of the test will be communicated by the testing laboratory to [INSERT LOCAL SYNLAB COUNTRY ENTITY] and to the physician mentioned on this form or its/their representatives.

Pursuant to the best practices and quality standards of clinical laboratories, the patient acknowledges that SYNLAB may use the leftover specimen and the patient's medical and genetic information, in an anonymized form (unless forbidden by applicable legislation) for research or quality assurance purposes. Such uses may result in the development of commercial products and services. The patient will not receive notice of any specific uses and will not receive any compensation for these uses. In any event, all such uses will be in compliance with applicable legislation.

Tick the box if you would not like your sample to be used for research purposes.

In compliance with Law 14/2007, of 3 July, on Biomedical Research, specifically its Articles 47 and 48, the prescribing physician must obtain informed consent to perform genetic tests. The patient's signature on this consent form is aimed at complying with this requirement.

Patient's name and surname(s): _____

Patient's or legal representative's signature: _____ Date: / / (day/month/year)

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60

BRCA+16GENES: Results

BRCA^{+16 GENES} results are reported as follows:

- ✓ **Pathogenic mutations:** Variants linked to disease (class V).
- ✓ **Probably pathogenic mutations:** Variants probably linked to disease (class IV).
- ✓ **Variants of uncertain significance:** Variants suspected of pathogenicity without decisive evidence (class III).
- ✓ **No mutations detected:** No mutations have been identified.

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BRCA+16 GENES

Patient report

Hereditary predisposition to breast, ovarian and endometrial cancer

PATIENT INFORMATION

First name:
Last names:
Gender:
Date of birth (Day/Month/Year):

ORDERING PHYSICIAN/CLINIC INFORMATION

Doctor:
Medical speciality:
Center/Hospital:

SAMPLE INFORMATION

Draw date:
Receipt date:
Sample ID:
Sample origin ID:
Report date:

REASON FOR THE STUDY

RESULT

NO MUTATION DETECTED

INTERPRETATION

Type here

No mutations classified as pathogenic, probably pathogenic or of uncertain significance have been identified in this study.

A result without pathogenic or probably pathogenic mutation detected does not eliminate the risk of breast, ovary and endometrial cancer.

This report is to be interpreted in the clinical context and family history of the patient.

BRCA+16GENES: Key points

- ✓ **BRCA+16 GENES** is a complete panel that **includes the genes contemplated by the NCCN guidelines in which specific patient management is established.**
- ✓ International experts recommend, based on experience, a **genetic screening for breast and ovarian cancer** that allows to stablish **specific prevention strategies for mutation carriers.** This type of screening is not contemplated in clinical guidelines for this purpose, but it is in certain patients with a history of risk.
- ✓ Aimed to the **healthcare providers** but **patients** can also be a target (high diffusion of the relationship between breast and ovarian cancer and genetic tests in the media).
- ✓ When a mutation is detected with **BRCA+16 GENES**, it is advisable to perform the analysis of the specific **mutation or deletion/duplication in the relatives.**
- ✓ **Genetic counselling** is key.

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- **Breast and ovarian cancer**
- **BRCA⁺¹⁶ GENES**
 - Genes panel and associated risk
 - Analysis technology
 - Recommendations from medical associations
 - Indications
 - Important considerations
 - Advantages
 - Sample collection kit
 - Test Requisition Form and Informed Consent
 - Results
 - Key points

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